

University of MN, Minnesota Nano Center

Standard Operating Procedure

Equipment Name: Oriel flood exposure system
Badger Name: oriel **Revision Number:** 3
Model: 8095 **Revisionist:** Paul Kimani
Location: Bay 4 **Date:** October 29, 2013

Description

The Oriel is a flood exposure system that produces a collimated and uniform beam of the UV light in the 350 - 450 nm region of the spectrum and is ideal for exposure of many photoresists. Most of the undesired light is filtered by internal optics.

1. Safety

- a. Avoid exposure to the direct, reflected or diffused UV light from the lamp.
- b. Wear UV glasses when the lamp is ON.
- c. Do not touch the lamp.

2. Restrictions/Requirements

- a. The volts/amperes cannot be adjusted.

3. Required Facilities

4. Definitions

5. Setup

- a. Check to see that the lamp is on by flipping the shutter toggle switch to **OPEN**.
- b. The lamp stays **ON**; the shutters for the lamp open and close exposing the sample to **UV** light. If the lamp is not on, press the black **START** button and wait **30 minutes** for the bulb to warm before continuing.
- c. The lamp current can be read either in Volts or Amperes. Move the switch to the desired method.

6. Operating Instructions

- a. Remove the black movable screen and set it aside.
- b. Load the sample, centering it on the middle of the X-sign.
- c. Replace the moveable black screen.
- d. Flip the shutter toggle switch to **OPEN** and start timing the exposure.
- e. Flip the shutter toggle switch to **CLOSE** when exposure time is complete.
- f. Remove the moveable black screen and remove the sample.
- g. Replace the moveable black screen.

a. Problems/Troubleshooting

- a. If the lamp does not turn on after restarting the lamp, it may need to be replaced.
- b. Notify the process staff.