

University of Minnesota Nano Fabrication Center

Standard Operating Procedure

Equipment Name: Varian E-beam evaporation

Coral Name:	ebevap-varian	Revision Number:	6
Model:	3118	Revisionist:	Kevin Roberts
Location:	Bay 3 - Keller Hall	Date:	2/10/2015

1 Description

The electron beam evaporator is used to evaporate solid dielectrics (granules, no powders) onto substrates. Evaporation is done under a high vacuum in a water cooled bell jar chamber. Evaporation is achieved by heating a source with an electron beam. As the source material evaporates, it forms a thin film on the samples.

2 Safety

- If you are evaporating and the building alarm sounds, TURN OFF THE POWER SUPPLY CIRCUIT BREAKER and leave immediately.
- To prevent sodium contamination, wear poly gloves whenever handling the source metals or the inside of the chamber.
- Wear UV glasses when viewing the beam to prevent eye damage.

3 Restrictions/Requirements

- Must be a qualified user
- Log in and out of the system using Badger
- Fill out the logbook.

4 Required Facilities

- Compressed air 60psi
- Process city water
- Exhaust

5 Definitions

- Ion gauge filament. Measures the pressure of the chamber while pumping with the diffusion pump (Cryo Pump)
- Hearth. Located inside the chamber and holds the metal sources.
- Planetarium. The fixture that holds the substrates inside the chamber.
- Shutter. A metal paddle that will cover/uncover the source metals.

6 Setup

Sample Loading

- Make sure Mode dial is on MAN.
- Turn off ion gauge using IG1 button.

University of Minnesota Nano Fabrication Center

Standard Operating Procedure

- c Switch HI VAC valve to CLOSE.
- d Switch AIR RELEASE under MANUAL VALVES label to OPEN.
- e Make certain "Convectron Gauge A" is above 7.0×10^{-2} Torr for 2 minutes.
- f Raise bell jar using HOIST RAISE switch.
- g Switch AIR RELEASE under MANUAL VALVES label to close.
- h Load samples and source(s) as you normally would.
- i Turn on the PROGRAM BOARD (the access panel needs to be flipped down). The T/X light should blink, then stop. If the light keeps blinking, the crystal needs to be replaced by a staff member.
- j Lower bell jar using LOWER HOIST switch. (Make sure black arrows are aligned).
- k Make sure AIR RELEASE is in closed position.
- l Switch ROUGH under MANUAL VALVES label to OPEN.
- m Watch Rough Pump Gauge ("A"). Allow it to reach 1×10^{-1} Torr (~10 min.).
- n Switch ROUGH to CLOSE position once 1×10^{-1} Torr is reached.
- o Switch HI VAC to OPEN position.
- p Wait 30 seconds after cross-over, then turn Ion Gauge on using "IG1" button.

7 Operating Instructions

- 1 The minimum pressure to operate is 6×10^{-5} Torr.
- 2 Enter the parameters on the program board. (this should already be on).
- 3 Make sure that the CONTROL POWER is set to AUTO.
YOU MUST PRESS ENTER after each parameter is entered.

University of Minnesota Nano Fabrication Center

Standard Operating Procedure

- 4 RISE TIME: is the time spent to ramp from zero power to SOAK POWER 1.
- 5 SOAK TIME 1: is the time spent at SOAK 1.
- 6 SOAK TIME 2: half of this time is spent ramping from SOAK POWER 1 to SOAK POWER 2, the other half is spent at SOAK POWER 2 .
- 7 RATE: this is the evaporation rate per second. Do not go above 5.
- 8 TOOLING A ratio for proper thickness measurement by the crystal. It is posted on the bell jar.
- 9 SOAK POWER 1: SOAK POWER 1 should be set 1 – 2% below SOAK POWER The purpose of both SOAK 1 and SOAK 2 is to evenly heat the source. This will allow the metal to out gas and also heat evenly, so the evaporation begins, it will be uniform.
- 10 SOAK POWER 2 should be set 1 – 2% above SOAK POWER 1.
- 11 MAX POWER: this setting prevents the beam from exceeding this level and damaging the equipment.
- 12 IDLE POWER: always set at zero (0).
- 13 GAIN: this controls the rate stability. The lower the gain, the more stable the deposition rate. Typical gains range from 1 to 10, (2) two is used most often.
- 14 DENSITY and ZRATIO are characteristics of the evaporant. See the process personal if you need to refer to the Film Evaporation Reference
- 15 SOURCE/ SENSOR: always set at 1/1.
- 16 Using blue rocker switch, Toggle out of Program Mode.
- 17 Wait for the bell jar to pump down to the desired vacuum 6×10^{-5} Torr.

POWER SUPPLY SET UP.

- 1 Make sure the crucible selector is set to the correct position.
- 2 Turn the circuit breaker to the POWER SUPPLY on.
- 3 Turn the SOURCE CONTROL POWER on.
- 4 All four interlocks should light GREEN. If not, the beam cannot be turned on: Water
Transformer
Vacuum
Doors (closed)
- 5 Turn the SWEEP CONTROLLER on, it is located on the back of the panel of the unit.
- 6 Turn the HIGH VOLTAGE on. Press the H.V. ON button.
- 7 Turn the FILAMENT CURRENT on. Press the Fill ON button.

University of Minnesota Nano Fabrication Center

Standard Operating Procedure

- 8 Set up the SWEEP PATTERN on the sweep controller
Patter should be set to 1 (circle)
Turn the sweep current by pressing ON.
Turn DC Bias ON.
Press PRM.

EVAPORATE

- 1 Turn the ROTATION on.
- 2 Turn Oxygen switch on (1 sccm)
- 3 Start the AUTO CYCLE. Press START, ZERO, START.
- 4 When the cycle reaches the DEPOSIT step, OPEN the SHUTTER.
- 5 After the desired thickness, when it has moved on to
Thickness 2 = 0, close the SHUTTER, and press STOP on the
Programmer.

SHUTDOWN

- 1 STOP will be lit. Turn off Programmer.
- 2 Press the H.V. RESET on the power supply
- 3 Press the FIL OFF on the Source Controller.
- 4 Turn the Circuit Breaker OFF on the POWER SUPPLY.
- 5 Turn the SOURCE CONTROLLER off.
- 6 Turn the SWEEP CONTROLLER off.
- 7 Turn the Oxygen off.
- 8 Turn the ROTATION off.
- 9 Turn the Ion Gauge off.
- 10 Allow source to cool for 5 minutes.
- 11 Switch HI VAC valve to CLOSE.
- 12 Switch AIR RELEASE under MANUAL VALVES label to OPEN.
- 13 Make certain "Convectron Gauge A" is above 7.0×10^2 Torr for 2 minutes.
- 14 Raise bell jar using HOIST RAISE switch.
- 15 Switch AIR RELEASE under MANUAL VALVES label to close.
- 16 Unload samples and source(s) as you normally would.

University of Minnesota Nano Fabrication Center

Standard Operating Procedure

- j Lower bell jar using LOWER HOIST switch. (Make sure black arrows are aligned).
- k Make sure AIR RELEASE is in closed position.
- l Switch ROUGH under MANUAL VALVES label to OPEN.
- m Watch Rough Pump Gauge ("A"). Allow it to reach 1×10^{-1} Torr (~10 min.).
- n Switch ROUGH to CLOSE position once 1×10^{-1} Torr is reached.
- o Switch HI VAC to OPEN position.
- p Wait 30 seconds after cross-over, then turn Ion Gauge on using "IG1" button.

8 Problems/Troubleshooting

- a The system is not pumping down. If the samples have moisture on them, i.e. water or photoresist, the chamber will take longer to reach the desired pressure.
- b If the samples are dry, there may be particles on the O-ring preventing a good seal. Vent the system and wipe off the O-ring with a wipe soaked with methanol and try to pump down the chamber again.
- c The T/X light keeps on flashing when the Program Board is turned on. This indicates that the crystal that monitors the evaporated metal needs to be changed by a staff member or the replaced crystal has been installed improperly.

University of Minnesota Nano Fabrication Center

Standard Operating Procedure

9 Special Notes on TiO₂

Depositing SiO₂ and Al₂O₃ vs Depositing TiO₂ — How They Differ

Al₂O₃ and SiO₂ work fine with crystal feedback, and dep. rate and final thickness can be controlled with crystal controller. Simply program the controller, and it will do everything for you.

TiO₂ does not work with the crystal controller, and the dep. rate will be erratic and uninformative. Also, final thickness displayed will be in error.

To deposit TiO₂ effectively:

- 1) Program TiO₂ with dep. rate 1-5 A/sec and give large (10.00 kA) final thickness to controller.
- 2) Program Power 1 = 5%, Power 2 = 10% and program Max. Power to either 10, 15, 20, or 25%
- 3) In practice the controller will go to this max. value and maintain that power while 'seeking' to reach input dep. rate of 1 to 5 A/sec.
- 4) Simply time run with amount of time shutter is open, and that will determine thickness based on the rates below.

We do not know why TiO₂ does not work with crystal control, but it is believed it may be that the crystal is simply too far away from the source.

TiO₂ Deposition Rates

With no oxygen flowing --

10% power = 27 A/min

20% power = 100 A/min

25% power = 140 A/min

With 1 sccm O₂ flowing, resulting in 9×10^{-5} Torr pressure --

25% power = 94 A/min

With 5 sccm O₂ flowing, resulting in 1×10^{-4} Torr pressure --

25% power = 75 A/min