

University of Minnesota Nano Center

Spinner CEE1- Standard Operating Procedure

Badger name:	K2 Spinner CEE1	Revision #:	2
Model:	200X	Revisionist:	Laura Parmeter
Location:	Keller- Bay 2	Date:	April 3, 2020

Table of contents

1. Introduction
2. Safety
3. Restriction/Requirements
4. Required Facilities
5. Definitions
6. Process Button
7. Operating Instructions
8. Problems/Troubleshooting

University of Minnesota Nano Center

Spinner CEE1- Standard Operating Procedure



1. Introduction

1.1. The CEE-1 spinner is a PC-controlled spinner with a touch screen interface and display used for applying photoresist uniformly on a substrate. It is capable of spin speeds from 0 to 12,000 rpm and spin speed accelerations of 0 to 30,000 rpm/s. Spin speeds and time can be altered to achieve desired photoresist thickness. Substrate size capability is from < 1cm to 200 mm round. Over 250,000 recipes with multiple steps can be programmed.

2. Safety

- 2.1.** For all chucks, ensure that the notch on the chuck is aligned to the drive pin on the spindle and push the chuck all the way down
- 2.2.** All chucks except the piece chuck need to be secured with the vented screw (with hole in center) that has been tightened with the provided torque hex wrench.

3. Restrictions/requirements

- 3.1.** Must be a qualified user on the CEE1 spinner
- 3.2.** Use only resist or material approved for use in spinner CEE-1. These are Shipley 1805, 1813, 1818, AZ 9260, SPR 220 7.0, **SPR 955 0.7CM**, Futurrex NR71 1500P, NR9 1500PY, PMMA in anisole and any other materials approved by NFC staff. All dispensed materials are held in one common waste storage tank
- 3.3.** Regularly used or **standard recipes are saved with 0_xxxx or 1_xxxx naming system.** These recipes will appear at the top of the list and no other recipe should precede them. Recipes that do not adhere to the above requirements will be deleted.
- 3.4.** Maximum spin speed (12,000 rpm) or maximum acceleration (30,000 rpm per sec) should not be exceeded

4. Required facilities

- 4.1.** Nitrogen supply (45psi – 55psi)

University of Minnesota Nano Center

Spinner CEE1- Standard Operating Procedure

- 4.2. Voltage range: single phase 100 – 120V AC, 10 amps
 4.3. Vacuum source: 25” Hg
 4.4. Exhaust: 50 cfm at 0.2 water



5. Definitions

- 5.1. **Recipe Name:** The name of the currently loaded recipe
 5.2. **Live update:** Not applicable to this spinner
 5.3. **Load:** Brings up the TODO Recipe Select Screen to select for processing
 5.4. **Edit Recipe:** This button will navigate to the Spin process Editor if a recipe is loaded. If a recipe is not loaded, it will go to the TODO Recipe Select Screen
 5.5. **Lid Closed:** A lamp showing the state of the Lid Closed sensor
 5.6. **Substrate Present/missing Indicator:** A lamp showing the state of the Spin Chuck Vacuum sensor
 5.7. **Step Indicator:** shows the current step of the process
 5.8. **Time Indicator:** Shows time remaining on the current step
 5.9. **Speed Indicator:** Shows the current speed of the spindle
 5.10. **Exhaust:** Shows the current settings of the programmable exhaust
 5.11. **Dispense:** Not applicable
 5.12. **Iterations:** Not applicable
 5.13. **Start Centering:** See below under **Process Button**
 5.14. **Vacuum hold/release:** this button allows the user to actuate the vacuum of the spin-chuck as they center the substrate.

University of Minnesota Nano Center

Spinner CEE1- Standard Operating Procedure

5.15. **Center:** Repeatabe centering. Use this button to check centering at any time.

6. **Process Button:** This has four modes:

6.1. **Start Centering:** Starts the wafer spinning very slowly to check for centering on the spin chuck

6.2. **Start Process:** Starts the selected process

6.3. **Abort:** Aborts a currently running process

6.4. **OK:** Turns off the process complete buzzer


6.5. **Vacuum Hold/Release button:** Allow users to actuate the vacuum of the spin-chuck as they center the substrate.

7. **Operating instructions**

7.1. Log into Badger and enable spinner CEE-1

7.2. On the main screen, tap the **Login** button (if previous user had logged out). This will bring up the keyboard to enter your password. User password is **1234**

7.3. Use appropriate size chuck. Available chucks are 5/16" for small samples, 4" for 100 mm wafers and an additional chuck for 150mm or 200mm wafers.




7.4. Press the **Run Spin Process** button, to bring up the spin process window

7.5. To run a spin process, load your substrate on an appropriate size chuck (substrate should fully cover all chuck surface), then load a recipe by pressing the **Load** button. *NB: If the*

University of Minnesota Nano Center

Spinner CEE1- Standard Operating Procedure

waste drain bottle is full, you will not be able to run any process until the waste bottle has been emptied.



- 7.6. A recipe **Selection Screen** will be brought up. Select a recipe by double clicking or pressing enter after your selection
- 7.7. Load wafer using provided centering tools (4" substrate loading depicted in the images below). Small substrates can be positioned without the use of centering tools. Press **Vacuum Hold** to actuate vacuum on the substrate.



- 7.8. Press **Start Centering**: The spin chuck spins for the 5 seconds to test for centering. Centering can be restarted by pressing the Center button. *NB: Centering happens whether the lid is open or closed.*
- 7.9. Dispense the desired amount of photoresist onto the center of the substrate. Close the lid.
- 7.10. Press **Start Process** to begin the spin-develop process. For the duration of the process the lid must remain closed. Once the process is started, the **Abort** button can be used at any time to halt the process.
- 7.11. **Process Complete**: Once the process is completed, an audible alarm may sound. Pressing the **OK** button will silence the alarm
- 7.12. **Clean the spinner bowl lining** using acetone and cleanroom wipes until there is no residual photoresist. Excess solvent will drain to the waste bottle via a drain hole below the chuck. Do not spray acetone on the chuck center where the screw is as this may ruin the motor.

University of Minnesota Nano Center

Spinner CEE1- Standard Operating Procedure

8. Problems/troubleshooting

8.1. Substrate present light not on: Ensure that CEE-1 is enabled on badger

8.2. Spin vacuum error: potential problems include vacuum lost during spin or spindle vacuum seal is worn out. Notify maintenance personnel.